

**LUDLUM MODEL 3001
MULTI-DETECTOR SURVEY METER
USER'S MANUAL**

November 2018

Serial Number: 25009185 and Succeeding

Firmware: n33 and Higher

LUDLUM MODEL 3001 MULTI-DETECTOR SURVEY METER USER'S MANUAL

November 2018

Serial Number: 25009185 and Succeeding

Firmware: n33 and Higher



LUDLUM MEASUREMENTS, INC
501 OAK STREET, P.O. BOX 810
SWEETWATER, TEXAS 79556
325-235-5494, FAX: 325-235-4672

STATEMENT OF WARRANTY

Ludlum Measurements, Inc. warrants the products covered in this manual to be free of defects due to workmanship, material, and design for a period of twelve months from the date of delivery. The calibration of a product is warranted to be within its specified accuracy limits at the time of shipment. In the event of instrument failure, notify Ludlum Measurements to determine if repair, recalibration, or replacement is required.

This warranty excludes the replacement of photomultiplier tubes, G-M and proportional tubes, and scintillation crystals which are broken due to excessive physical abuse or used for purposes other than intended.

There are no warranties, express or implied, including without limitation any implied warranty of merchantability or fitness, which extend beyond the description of the face there of. If the product does not perform as warranted herein, purchaser's sole remedy shall be repair or replacement, at the option of Ludlum Measurements. In no event will Ludlum Measurements be liable for damages, lost revenue, lost wages, or any other incidental or consequential damages, arising from the purchase, use, or inability to use product.

RETURN OF GOODS TO MANUFACTURER

If equipment needs to be returned to Ludlum Measurements, Inc. for repair or calibration, please send to the address below. All shipments should include documentation containing return shipping address, customer name, telephone number, description of service requested, and all other necessary information. Your cooperation will expedite the return of your equipment.

**LUDLUM MEASUREMENTS, INC.
ATTN: REPAIR DEPARTMENT
501 OAK STREET
SWEETWATER, TX 79556**

**800-622-0828 325-235-5494
FAX 325-235-4672**

LUDLUM MEASUREMENTS, INC.

501 OAK ST. / P.O. BOX 810

SWEETWATER, TEXAS 79556

Phone: 325-235-5494 800-622-0828 (USA)

Fax: 325-235-4672 E-Mail: ludlum@ludlums.com

Website: <http://www.ludlums.com> FEIN 75-1085764



DESIGNER AND MANUFACTURER

OF

*Scientific and Industrial
Instruments*


Model 3000 and 3001 Firmware Differences

October 2018


The two-alarm system mentioned in Section 4 of the Model 3000 or Model 3001 manual only applies to instruments with Firmware Version n33 and higher. For firmware below Version n33, there is only one alarm.

REVISION HISTORY

REV	DESCRIPTION	DATE	BY
1	VALID	5/26/13	DDW
2	ADD NOTES	1/25/14	DDW
3	ADD CABLE CLAMP	4/16/14	DDW
4	CHANGE COAX ROUTING	12/3/14	DDW
5	ECF# 3778	10/7/15	WJM
6	added 20-9008 EC3951	1/8/16	JGW
7	CHANGED HANDLE ECF 4925	10/25/17	JWI



DETAIL Y
SCALE 1 : 1


DWN	DATE	CHK	DATE	APP	DATE
JWI	10/25/17			JWI	11-13-17
DWG NUM: 4498-409				SCALE: 1:1	
TITLE M 3001 FRONT PANEL					
 LUDLUM MEASUREMENTS, INC. 501 OAK STREET SWEETWATER, TEXAS 79556			SERIES 498	SHEET 409M	

JWI

Table of Contents

<i>Introduction</i>	1
<i>Getting Started</i>	2
Unpacking and Repacking	2-1
Battery Installation	2-1
Instrument Operational Test	2-2
Sigma Audio	2-4
Detector Failure Diagnostic	2-5
Detector Over Range	2-5
Detector Overload	2-6
Instrument Use and Controls	2-7
RATE Mode Operation	2-9
MAX Mode Operation	2-10
COUNT Mode Operation	2-11
<i>Specifications</i>	3
<i>Setup Mode</i>	4
Setup Overview	4-1
Setup Parameters	4-2
Entering Setup Mode	4-4
Setup Mode Operation	4-5
<i>Datalogging</i>	5
Datalogging Operation – Mode 1	5-2
Datalogging Operation – Mode 2	5-3
Datalogging Operation – Mode 3	5-4
Specifications	5-4
<i>Software</i>	6
Connecting to Software	6-1
Extended Features	6-2
<i>Advanced Features</i>	7
Dead Time Correction	7-1

Units	7-2
Software Calibration Tools	7-3
Other Device Data	7-3
Real-Time Streaming	7-4
<i>Safety Considerations</i>	<i>8</i>
Environmental Conditions for Normal Use	8-1
Warning Markings and Symbols	8-1
Cleaning and Maintenance Precautions	8-2
<i>Revision History</i>	<i>9</i>
<i>Recycling</i>	<i>10</i>
<i>Options</i>	<i>11</i>
<i>Standard Parts List</i>	<i>12</i>
<i>Auxiliary Communications</i>	<i>A</i>
AuxComm Overview	A-1
Requirements	A-1
Usage	A-2
Settings	A-5

A dark gray square graphic with the word "Section" in white at the top and a large white number "1" in the center.

Introduction

The Model 3001 is an ergonomic, lightweight instrument which can be used for radiation survey purposes. It features the ability to measure radiation in count rate, exposure rate/dose, activity rate, integrated exposure/dose, time-averaged rates, and scaler counts. The Model 3001 supports up to four external radiation detectors to detect alpha, beta, or gamma radiation. Each detector may have its own set of calibration and user parameters. The DETECTOR button is used to select the active detector.

The instrument features a large backlit LCD (liquid crystal display), an audio warning, and easy, intuitive use. The unit body is made of lightweight but durable plastic. It is intended for outdoor use and can resist splashing water. The foam grip provides a comfortable, ergonomic grip but may be removed for ease of decontamination.

Three modes of operation are available for the Model 3001 – RATE, MAX, and COUNT. RATE mode operation will display the current count, exposure, or activity rate. MAX mode is used to capture the highest count, activity, or exposure rate detected – useful for finding a peak rate, or frisking when the display is not visible. Two sets of units (primary and secondary) for RATE and MAX modes can be chosen from among cps, cpm, Bq, dpm, R/h, or Sv/h. The user can switch between these two units by simply pressing the UNITS button.

COUNT mode allows the user to perform a count for a predetermined time. Depending on the count units chosen, the result can be a scaler count (in counts or disintegrations), a time-averaged rate (cpm, dpm, Bq, cps), a time-averaged exposure or dose (R/h, Sv/h), or an integrated exposure or dose (R or Sv).

The standard handle includes the LOG pushbutton that allows the instrument to log up to 1000 data points when activated with the optional Lumic Datalogging Kit.

The display will be automatically backlit if light levels are low. The display backlight can also be configured for "Continuous On" operation. RATE, MAX, and COUNT modes can be silent or utilize a "click" audio; the "click" audio is always silent during COUNT mode by default. A "sigma" audio mode can be enabled (disabling the "click" audio) in RATE and MAX modes when scintillator detectors are used, and this mode makes it easy for the user to find small increases above the background radiation level. In this mode, the instrument measures background for 8 seconds after power-up, and then beeps whenever the rate increases by a small amount. Users are freed from watching the numeric value to "find" something. They can simply listen for multiple beeps.


Setup of the instrument is accomplished through the front-panel buttons, or through software available from Ludlum Measurements.

Front-panel setup can be disabled via the internal switch on the Model 3001 in order to protect settings from inadvertent changes. When the settings are protected, the user can only view the settings, and the password-protected software must be used to change any of the settings. The Model 3001 can be configured through software to allow changes to the alarms and or count time even when in protect mode.

The unit is operated with four alkaline AA batteries for operation from -20 to 50 °C (-5 to 122 °F). Battery life is approximately 500 hours under normal usage. A low-battery indicator on the LCD warns when less than 16 hours of battery life remain.

The Model 3001 may optionally be provided with a 1/8-inch stereo headphone jack that would allow for headphone use for audio output. It is located to the left of the detector connector.

For more advanced use: See Appendix A for information on our Auxiliary Communications features, which allow the instrument to expand its capabilities with a variety of external devices through a standard serial interface.

A dark gray square with the word "Section" in white at the top and a large white number "2" in the center.

Getting Started

Unpacking and Repacking

Remove the calibration certificate and place it in a secure location. Remove the instrument and ensure that all of the items listed on the packing list are in the carton. Check individual item serial numbers and ensure calibration certificates match between instruments and detectors (if applicable). The Model 3001 serial number is located on a label on the front side of the unit.

To return an instrument for repair or calibration, provide sufficient packing material to prevent damage during shipment.

Every returned instrument must be accompanied by an **Instrument Return Form**, which can be downloaded from the Ludlum website at www.ludlums.com. Find the form by clicking the "Support" tab and selecting "Repair and Calibration" from the drop-down menu. Then choose the appropriate Repair and Calibration division where you will find a link to the form.

Battery Installation

A low-battery indicator appears at the top of the LCD when less than 16 hours of battery life remain. When this indicator is present, follow these steps to replace the four standard AA batteries:

1. Turn the instrument over so that the bottom of the instrument is facing up.
2. Use a straight, medium-sized screwdriver to turn the single screw on the battery cover one-quarter turn counter-clockwise.
3. Release and remove the battery cover.
4. Replace all four AA batteries.
5. Replace the cover and turn screw one quarter of a turn clockwise to secure.

Instrument Operational Test

Turn the instrument ON by pressing the ON/ACK button for about a second, and then releasing.

The instrument should activate all the LCD segments and the audio. Observe the device during this time. If any LCD segments are missing, or audio fails to work, the device is in need of repair. Please refer to Figure 1 below.




Figure 1: Startup display for Model 3001, with all LCD segments shown.

The instrument then displays the firmware version. Please refer to Figure 2 below.



Figure 2: Firmware version display.

Note: Firmware versions with the Bq units are equipped with auxiliary communication.

The instrument then displays the battery voltage. Please refer to Figure 3.




Figure 3: Battery Voltage display

The instrument then displays the number of stored records if datalogging is enabled. Please refer to Figure 4 below.



Figure 4: Startup display for Optional Data Logging firmware, showing 155 stored records.

The instrument will then move to normal operation, displaying the current rate for the Primary units (default: mR/hr).

If the Sigma Audio option is selected, the unit will display a countdown from :08 to :01 (in seconds) as the unit measures background radiation levels.

The user may select the Secondary units (default: cpm) by pressing the Units button.

Ensure that the low-battery indicator is not present. If the low-battery indicator is present, replace the batteries as soon as possible. Should the instrument detect a battery voltage that is high enough to power on, but too low to safely operate, the display will blank and the low-battery icon will flash. Normal operation will not be available until the batteries have been replaced. Under extreme low-battery conditions, be aware that the unit may not even turn on or may turn itself off abruptly.

A reference reading with a check source, 1 μCi (37 kBq) of ^{137}Cs for example, should be obtained at the time the instrument is received in the field. Small check sources of radiation are available from Ludlum Measurements. While exempt from many regulations because of their small size, these sources are large enough to produce a response on this instrument. If this procedure is done routinely with the same radiation source, instrument malfunction may be detected when anomalous readings are observed. If at any time the instrument fails to read within 20% of the reference reading when using the same check source, it should be sent to a calibration facility for recalibration and/or repair.

Example log reading:

Check Source # _____ Rate _____ Units _____

Once this procedure has been completed, the instrument is ready for use.

Sigma Audio

The Ludlum Model 3001 has the standard "click" audio that is widely used by radiation instruments. However, it also has a Sigma Audio mode that can be enabled to change the audio sound. This mode is useful for scintillation detectors because their high count rate makes the "click" audio less useful.

In the sigma audio mode, the instrument measures the background radiation alarm for eight seconds on power-up, and then automatically sets a low threshold rate above the background rate. Then the instrument will produce a beep any time that the instrument "sees" more than this threshold value. The user doesn't need to watch the display to find a small increase over background; the user may just listen for multiple beeps. Thus, the unit has an audio alarm that, on power-up, adjusts to just above the current background level and provides a sensitive audio indication to the user.

Note that this audio alarm can also work in conjunction with the fixed alarm, i.e. the user can have both a floating audio alarm (resulting in audio beeps) based on the background level, and a fixed tone audio and a steady ALARM icon when a predetermined fixed alarm level is exceeded.

Detector Failure Diagnostic

Note that the Model 3001 has its own diagnostic tests to ensure that the detector is functioning correctly. The Model 3001 can detect when the radiation detector is malfunctioning and will flash the display to indicate a fault. If the detector stops detecting radiation for a settable number of seconds, the Model 3001 will flash a zero reading for the currently selected units. This indication is common if the unit is powered up without a detector connected. If this indication is observed with a connected detector, remove the unit from service and have it evaluated by a qualified repair and calibration technician.




Figure 5: Detector Failure display (shown for cpm); will also flash.

Detector Over Range

If the detector is exposed to high levels of radiation or has an internal malfunction that causes it to count high or excessively, the unit flashes the maximum rate (999) for the currently selected units as a warning. The icon ALARM is also displayed at the top of the screen. The user should ensure whether this condition is being caused by a high radiation field, by a shorted cable, or by internal malfunction. With some detectors, this display could be an indicator of a light leak, possibly caused by a puncture or tear in the detector face. If a Maximum Display parameter is set, then it is possible for the display to flash the Maximum Display value under this condition.




Figure 6: Detector Over Range (shown for kcpm); will also flash.

Detector Overload

As another diagnostic test, the Model 3001 monitors the HV supply's detector current. A current overload threshold can be set via Setup Mode. (A setting of 0 disables this alarm.) In general, this alarm setting can be used to detect when there is some detector failure, a cable failure, or a possible high level of radiation. When this alarm is triggered, the LCD will display OL and an alarm tone will sound. Once the detector current goes below the threshold, the Model 3001 will return to normal operation.



Figure 7: Detector Overload Alarm; will also flash.

Instrument Use and Controls

With four front-panel buttons (as seen in front-panel drawing in front of manual), the Ludlum Model 3001 is simple and easy to use with minimal training required. Default operation is RATE mode, and the display shows the current count rate using the Primary units. Long-pressing the UNITS button will switch between Primary and Secondary units. Long-pressing the MODE button will switch the instrument to MAX mode, which will display the highest count rate detected. Long-pressing the MODE button again will switch it to COUNT mode, which will display the COUNT timer. Note that either or both of the MAX and COUNT modes can be locked out in the setup process.

ON/ACK button: Used to power the Model 3001 ON and OFF, reset MAX mode, start/reset the COUNT Timer, and acknowledge audio alarms.

- Power On: Press for approximately one second and release (all LCD segments will activate, and firmware version will be shown).
- Power Off: Press and hold for approximately four seconds. The display will show a 3, 2, 1 countdown for the final three seconds of shutdown. Releasing the ON/ACK button during shutdown will return the device to the previous state of operation. At completion of the shutdown count, the LCD will go blank.

Short Tap (tap-pressing the button and releasing in under a half second):

- Normal Operation: Will acknowledge/silence alarms in all modes of operation.
- Change Log Location Selected Position (if data logging is enabled).
- Will move selection position in the device menu.

Long Press (pressing the button and releasing for at least a half second):

- Start/reset COUNT Timer in COUNT mode.
- When Sigma Audio is enabled, pressing the ON/ACK button will take a new background reading and update the Sigma Audio level (in Rate Mode only).
- Reset/Zero: When in Rate mode, a long press of the ON/ACK. Will reset/zero the averaged rate: the minimum display will be shown before returning to the exposure rate. When Sigma Audio is enabled, Rate Reset is disabled.
- Will reset MAX mode display.

UNITS button:

Short Tap (pressing the button and releasing under a half second)

- Used to switch between Primary and Secondary units in all available modes.

DETECTOR button:

Short Tap (pressing the button and releasing under a half second)

- Used to select the active detector, as evidenced by the corresponding DET1, DET2, DET3, or DET4 LED. Unit saves this setting prior to shutting down, so it will power up with the same detector setting.

Long Press (pressing the button and releasing for at least a half second)

- Only used in the device menu. While viewing the Detector Overload parameter, a long press will allow the user to view the actual current level. Another long press will return user to overload setting.

MODE button:

Short Tap (pressing the button and releasing under a half second)

- Used to advance between the three operating modes: RATE, MAX, and COUNT. Note that MAX and/or COUNT mode may be disabled from use.
- The device menu increments selected digit/text to the next available option.

RATE Mode Operation

In RATE mode, the current count rate will be displayed.

Long-pressing the ON/ACK button for at least a half second will reset the averaged rate before releasing the button. The Reset feature must be enabled through Lumic Calibration software (see parameter: RateResetBtn) and will be disabled when Sigma Audio is enabled.

Tapping the UNITS button will switch the displayed value between the Primary and Secondary Units.

Tapping ON/ACK button will turn the “click” audio on/off. When Sigma Audio is enabled, tapping the ON/ACK button will toggle the Sigma Audio beep on/off. If dual Level audio is enabled, it will toggle between “Hi, LO, and OFF.”

If an alarm condition is present, tapping the ON/ACK button will acknowledge and turn off the continuous tone alarm audio. Under an alarm condition, the ALARM display indicator will remain on. Alarms are non-latching in RATE mode. When Sigma Audio is enabled, and an alarm condition is not present, a long press (at least a half second) of the ON/ACK button will reset the Sigma Audio alarm level. The Sigma count can only be reset in RATE mode.




Figure 8: RATE mode display showing typical background radiation rate and the low-battery icon.

MAX Mode Operation

While in MAX mode, the highest detected count rate (since the last reset) is displayed. The word MAX will be displayed when in MAX mode.

Tapping the UNITS button will switch the displayed value between the Primary and Secondary Units.

Tapping the ON/ACK button will turn the “click” audio on/off. When Sigma Audio is enabled, tapping the ON/ACK button will toggle the Sigma Audio beep on/off.

Under a non-alarm or alarm condition, a long press of the ON/ACK button will reset the display. When Sigma Audio is enabled, the Sigma Audio alarm level cannot be reset in Max Mode.

If an alarm condition is present, tapping the ON/ACK button once will acknowledge and turn off the continuous tone alarm audio. (The “click” audio will remain as selected under non-alarm conditions.) Tapping the ON/ACK button a second time will toggle the audio setting. Under an alarm condition, the ALARM display indicator will remain on until the display is reset. Alarms in MAX mode latch with the display.

If other operational modes are available, tapping the MODE button will move to the next available operational mode.




Figure 9: MAX mode operation display with ALARM indicator.

COUNT Mode Operation

When entering COUNT Mode from another operational mode, the currently selected COUNT Unit will be displayed for approximately one second. The purpose of COUNT mode is to count for a predetermined amount of time, and to display the results on the display. Note that the predetermined count time can be from 1 second to 10 minutes, or can be set to zero to enable continuous counting until stopped by the user.

Count mode operation is very flexible, depending on the units chosen. A common choice is for the count mode to just perform a scaler count for a specified time, with a resulting answer in counts (equaling detected radiation events).

There are two different options for COUNT mode, which can be set through software.

Go to “Device Count DisplayMode” in the Device or Others tab (depending on software version). Enter the desired option into the “Value” field.

- 0 - will only show the countdown timer while the timer is active (default selection).
- 1 - will cycle between showing the countdown timer and the current reading.
- 2 - will only show the current gathered reading.

Note:

Tapping the LOG button will start a count and automatically log the result. A long press of the LOG button during the count will toggle the display between the timer and reading. Audio ‘clicks’ are disabled in COUNT mode by default

Device Count Mode AudioMode will enable\disable click audio during a COUNT

- 0 – Disabled (Mute icon will always be illuminated in COUNT mode.)
- 1- Enabled

Note:

Sigma audio is disabled in Count Mode. If Sigma Audio and Count Mode Audio are enabled, click audio will be utilized in Count Mode. If a result in terms of activity is desired, the scaler count can also be in units of "d" or disintegrations. But if the count units are chosen to be cpm or cps, then the resulting answer is an averaged count rate over the time interval. Similarly, if count units of Bq or dpm are chosen, the resulting answer is an averaged disintegration rate.

Note:

If the user desires the instrument to show results in terms of disintegrations/area (eg. dpm/100cm² or Bq/cm²), then the appropriate factor should be placed in the Efficiency parameter.

Other choices are to have count mode units of R/h or Sv/h, in which case the COUNT mode result is an averaged exposure or dose rate. But if count mode units of R or Sv are chosen, the result is shown in accumulated exposure or accumulated dose over the chosen count time. The following tables lists the possibilities:

UNITS	RESULT
c	counts per count time
d	disintegrations per count time
cpm, cps	count rate averaged over the count time
dpm, Bq	disintegration rate, averaged over the count time
R/h, Sv/h	exposure or dose rate, averaged over the count time
R, Sv	integrated exposure or dose over the count time

In COUNT mode, operation depends on the current state of the Count Timer.

When the Count Timer is Ready:

- The display will show the Count Time.
- Tapping the UNITS button will switch between the Primary and Secondary Count Units. The newly selected Count Units will be displayed for approximately one second, and the display will then return to the Count Timer.
- Long-pressing the ON/ACK button starts the Count Timer.
- If other operational modes are available, tapping the MODE button will move to the next available operational mode.

When the Count Timer is active:

- The display will show either Count Time (default), both Count Time and Counts, or just counts, depending on Count Display Mode.
- Doing a long press (pressing the button and releasing after at least a half second) the ON/ACK button will reset the Count Timer.
- Tapping the UNITS button will display the other unit if enabled.
- Click audio, if enabled, can be turned on and off by tapping the ON/ACK button.
- If an alarm condition occurs, the ALARM display indicator will turn on. Alarm audio will sound after the count is completed if the result is over the alarm threshold. (If the count unit is an averaged rate, it is possible for the ALARM icon to come on during the count, but if the averaged rate falls below the alarm threshold by the end of the count, no ALARM will be active.)
- If other operational modes are available, tapping the MODE button will cancel the current Count Timer and move to the next available operational mode.

When the Count Timer has finished:

- The display will show either the accumulated total for c, d, R, and Sv, or the timed ratemeter average for cps, cpm, Bq, dpm, or the average exposure or average dose in R/h and Sv/h.
- Tapping the UNITS button will switch between the Primary and Secondary Count Units.
- Pressing the ON/ACK will reset the Count Timer.

- If an alarm condition occurred during the Timed Count, a continuous audio tone will sound, and the ALARM display indicator will already be on. Tapping the ON/ACK button once will acknowledge and turn off the continuous tone alarm audio. Doing a long press (pressing the button and releasing after at least a half second) to the ON/ACK button will clear the alarm condition and reset the Count Timer. Alarm audio will sound after the count is completed if the result is over the alarm threshold.
- If other operational modes are available, tapping the MODE button will move to the next available operational mode.



Figure 10: COUNT mode operation showing COUNT Timer of 5 minutes, 30 seconds.

Section

3

Specifications

Detector: may be Geiger-Mueller (GM), scintillator, or proportional

HV Range: 400-1500 Vdc

Threshold Range: 4-100 mVdc

Resolving Time: approximately 5 μ sec as defined by IEC 60325

Alarms: count rate, exposure/dose, and count alarm setpoints adjustable over the display range

Overload: high count rate saturation protection prevents false display of lower count rates

Zero Protection: after a user-settable number of seconds of no pulses from detector, unit will flash a zero reading and the alarm audio will be triggered

LCD Display: 3-digit LCD with large 20 mm (0.8 in.) digits, (k)cps, (k)cpm, (k)Bq, (k)dpm, (μ)(m)R/h, (μ)(m)Sv/h, low-battery indicator, MAX, ALARM, AUDIO

Display Range:

- 0.0 cps to 99.9 kcps
- 0.00 cpm to 999 kcpm
- 0.00 Bq to 99.9 kBq
- 0.00 dpm to 999 kdpm
- 0.00 μ R/h to 999 R/h
- 0.00 μ Sv/h to 999 Sv/h

Headphone Jack (optional): 1/8 inch stereo jack for use with headphones for audio output. It is located to the left of the detector connector.

Backlight: built-in ambient light sensor automatically activates low-power LED backlight, unless internal dip switch is set to continuous-on (will reduce battery life)

User Controls:

- ON/ ACK – long press to turn ON, tap to acknowledge alarms and silence alarm tone, long press to reset Sigma Audio alarm, long press to reset Rate, and hold for OFF. Tap to turn “click” audio on/off, or turn Sigma Audio beep on/off
- MODE – alternates between RATE (count rate), MAX (captures peak rate), and COUNT (user-selectable preset count time from 0 to 10 minutes)
- DETECTOR – tap to change active detector to the next enabled detector
- UNITS – changes the units from count rate (cpm, cps), to dose/exposure ($\mu\text{Sv/h}$, mR/h) or disintegration (dpm, Bq)

Response Time: user-selectable from 1 to 60 seconds, or auto-response rate FAST or SLOW

Audio: click audio greater than 65 dB at 0.6 m (2 ft), alarm audio at greater than 72 dB, approximately 4.5 kHz

Power: four alkaline “AA” batteries

Battery Life: approximately 500 hours of operation (as low as 100 hours with backlight configured for continuous-on), 16-hour low-battery warning

Maximum Current: 35 mAdc

Construction: high-impact plastic with water-resistant rubber seals and separate battery compartment

Temperature Range: -20 to 50 °C (-5 to 122 °F), may be certified for operation from -40 to 65 °C (-40 to 150 °F)

Environmental Rating: NEMA (National Electrical Manufacturers Association) rating of 4x or IP (Ingress Protection) rating of 53

Size: 16.5 x 11.4 x 21.6 cm (6.5 x 4.5 x 8.5 in.) (H x W x L)

Weight: 0.9 kg (2.0 lb)

Section

4

Setup Mode

Warning!

Only advanced users or administrators should consider changing any of the parameters in the following section. Incorrect settings could jeopardize the safety of users depending on this instrument.

Setup Overview

Your instrument has been shipped from Ludlum Measurements only after passing electronic checkout, a 24-hour burn-in process, and a careful calibration process. Calibration papers are supplied with each instrument shipped from Ludlum Measurements.

Recalibration should be accomplished after maintenance or adjustments have been performed on the instrument. Recalibration is not normally required following instrument cleaning or battery replacement. Recalibration does not require any special tools or software to perform.

Ludlum Measurements offers a full-service repair and calibration department. Not only do we repair and calibrate our own instruments, we also service most other manufacturers' instruments. Calibration procedures are available upon request for customers who choose to calibrate their own instruments.

Note:

Ludlum Measurements, Inc. recommends recalibration at intervals no greater than one year, assuming that regular operational checks are performed. Check the appropriate local, state, and federal regulations to determine required recalibration intervals.

Setup Parameters

Setting	Setup Parameter	Model 133-6	Model 44-2	Alpha or Beta	Model 42-31H	Model 44-132
P1-1	Calibration Constant Mantissa	122	100	001	600	940
P1-2	Calibration Constant Exponent	04	08	01	04	8
P1-3	Dead Time Correction 1	36	10	5	5	5
P1-4	Dead Time Correction 2 Mantissa	360	TBD	TBD	TBD	TBD
P1-5	Dead Time Correction 2 Exponent	-12	TBD	TBD	TBD	TBD
P1-6	Efficiency	100	15	Alpha: 20 Beta: 15	15	15
P1-7	High Voltage Setpoint	550	800	800	1200	800
P1-8	Discriminator Threshold	35	35	35	2	10
P1-9	Detector Current Overload Threshold	0 (disabled)				
P1-10	Loss of Count Alarm Time	999	60	Alpha:600 Beta:60	0 (disabled)	60
P1-11	Detector Enable	En				
P2-1	Unit 1 - Units and Minimum Display	00.0 mR/h	000 μ R/h	000 cpm	0.00 μ Sv/h	00.0 μ R/h
P2-2	Unit 1 - Units Maximum Display	999 R/h	50.0 mR/h	999 kcpm	100.0 mSv/h	2.50 mR/h
P2-3	Unit 1 - RATE/ MAX Mode Alarm 1	0 (disabled)				
P2-4	Unit 1 - RATE/ MAX Mode Alarm 2	50.0 mR/h	5.00 mR/h	5.00 kcpm	50.00 μ Sv/h	500 μ R/h
P2-5	Unit 1 - Count Units and Minimum Display	00.0 mR/h	0.00 mR	000 counts	0.00 μ Sv	00.0 μ R/h
P2-6	Unit 1 - Count Units Alarm 1	0 (disabled)				
P2-7	Unit 1 - Count Units Alarm 2	000 R (disabled)	000 R (disabled)	0 counts (disabled)	000 Sv	000 R (disabled)
P3-1	Unit 2 - Units and Minimum Display	000 cpm	000 cpm	oFF	000 cpm	000 cpm
P3-2	Unit 2 - Units Maximum Display	999 kcpm	999 kcpm		999 kcpm	999 kcpm
P3-3	Unit 2 - RATE/ MAX Mode Alarm 1	0 (disabled)				

Setting	Setup Parameter	Model 133-6	Model 44-2	Alpha or Beta	Model 42-31H	Model 44-132
P3-4	Unit 2 - RATE/ MAX Mode Alarm 2	5.00 kcpm	5.00 kcpm	5.00 kcpm	5.00 kcpm	5.00 kcpm
P3-5	Unit 2 - Count Units and Minimum Display	000 counts	000 counts	oFF	000 counts	000 counts
P3-6	Unit 2 - Count Units Alarm 1	0 (disabled)				
P3-7	Unit 2 - Count Units Alarm 2	0 counts (disabled)	0 counts (disabled)		0 counts (disabled)	0 counts (disabled)
P4-1/2	Response Time/Auto Response Rate	S				
P4-3	Enabled Modes	Rate/Max/Count				
P4-4	Count Time	60				
P4-5	Auto Shutdown Time	0 hours (disabled)				
P4-6	Sigma Audio Mode	oFF				on
P4-7	Dual Level Audio	oFF				
P5-1	Datalogging Mode	1				
P5-2	Real Time Clock Month	1				
P5-3	Real Time Clock Day	7				
P5-4	Real Time Clock Year	14				
P5-5	Real Time Clock Hour	16				
P5-6	Real Time Clock Minute	30				
P5-7	Real Time Clock Second	00				

Entering Setup Mode

To enter Setup Mode, power down the Model 3001, then turn the unit back ON. Following the display of the Firmware version, when the instrument has begun normal operation, tap the MODE button three times (within four seconds) to enter Setup mode.

Note:

This process is different if you are in Sigma Audio Mode, rather than in the “click” Audio Mode. While trying to enter Setup Mode with Sigma Mode enabled, once the screen displays the firmware number and then begins the Sigma countdown, press the Mode button three times. Do not wait until the countdown is complete, which will prevent you from getting into the Setup Mode. There is no countdown for the “click” Audio Mode.

Entry to Setup Mode can be confirmed when the numeric portion of the display shows P-1, indicating the first setup page is selected. If you simply wish to view the parameters, select the desired Setup Page by tapping the MODE button. Tap the UNITS button to advance through the parameters available on the selected Setup Page. To return to normal operation, advance back to the Setup Page selection by pressing the UNITS button for at least a half second. Pressing the UNITS button again while the page selection is active will return to the **Setup Detector Selection**. Pressing the UNITS for at least a half second again will exit, or holding the UNITS button for approximately 5 seconds will exit from anywhere in the menu.

Note:

If the menu is not exited properly, any changes made will not be saved.

SETUP PROTECT: The Model 3001 parameters can be protected from unauthorized changes via the internal switch located on the Model 3001 circuit board. To change the switch, open the battery compartment and remove the batteries from the Model 3001. Next, loosen the six pan head screws that fasten the bottom cover.

Gently remove the bottom cover of the instrument. The DIP (dual in-line position) switch should now be visible in the upper left-hand corner of the circuit board.

To protect the Model 3001 from changes in Setup mode, slide DIP Switch 2 (the rightmost switch) to the ON (forward) position. If DIP Switch 2 is in the OFF (back) position, changes are allowed in Setup mode. Once the DIP Switch is set as desired, gently replace the back cover and the six pan head screws. Install the batteries, and replace the battery cover.

There are three different options for protect mode, which can be set through Lumatic Calibration Software. Locate the Parameter "Device Setup Prtct" in the Others tab. Enter the desired option into the "Value" field.

- **Normal:** All Parameters are protected from change through the device menu.
- **Bypass Alarms:** All Parameters except for Rate and Count Alarms are protected from change through the device menu.
- **Bypass Alarms Count:** All Parameters except for Rate Alarms, Count Alarms, and count time are protected from change through the device menu.

Note that with the DIP Switch 2 in the ON position. Setup mode may be entered and parameters viewed, but changes cannot be made.

DISPLAY BACKLIGHT 'Continuous On': The Model 3001 display backlight can be set to remain on continuously during operation. Follow the steps above for **SETUP PROTECT**, but use DIP Switch 1 for display backlight selection. Setting DIP Switch 1 to the ON (forward) position will configure the display backlight to remain on during operation. Set DIP Switch 1 to the OFF (back) position, and the display will be backlit only when light levels are low.

Note:

Setting the display backlight for continuous-on operation can result in reduced battery life. The backlight can be configured to always be off as well through the device menu or Lumatic Calibration Software.

Setup Mode Operation

Once the Model 3000 is in Setup mode, the Setup Page selection will be displayed on the LCD, and the Page number will be blinking, indicating it as the

selected item. The number of available parameters per Setup Page will be displayed using the graph-the number of segments indicating the number of parameters. Tap the MODE button to choose the Setup Page. Once the desired Setup Page is shown, tap the UNITS button to move to the first parameter of that Setup Page. Long-pressing the UNITS button will exit the menu.



Figure 11: Setup Page Selection display (showing page 4).

Long-pressing the UNITS button, no matter what parameter is shown, will first return the user back to the Page Selection screen and then exit if held for approximately 5 seconds in total. Exiting can also be accomplished by one or two individual long presses (at least a half second) of the UNITS button depending on where the user is in the menu. (If on the page selection, only one press will exit, but if viewing an individual parameter, two presses will be required.).

The list below shows the four setup pages and the parameters, in order, on each page.

Model 3001 List of Parameters (in order)

Page 1 (P-1)

- Calibration Constant Mantissa
- Calibration Constant Exponent
- Dead Time Correction 1
- Dead Time Correction 2 Mantissa
- Dead Time Correction 2 Exponent
- Efficiency
- High Voltage Setting
- Pulse Threshold
- Detector Current Overload Threshold
- Loss of Count Alarm Time

Page 2 (P-2)

- Primary Units and Minimum Display
- Primary Units Maximum Display
- Primary Units RATE/MAX Mode Alarm 1
- Primary Units RATE/MAX Mode Alarm 2
- Primary COUNT Units and Minimum Display
- Primary COUNT Units Alarm 1
- Primary COUNT Units Alarm 2

Page 3 (P-3)

- Secondary Units and Minimum Display
- Secondary Units Maximum Display
- Secondary Units RATE/MAX Mode Alarm 1
- Secondary Units RATE/MAX Mode Alarm 2
- Secondary COUNT Units and Minimum Display
- Secondary COUNT Units Alarm 1
- Secondary COUNT Units Alarm 2

Page 4 (P-4)

- Response Time
- Auto Response Rate
- Operational Modes
- Count Time
- Auto Shutdown Time
- Backlight Threshold
- Sigma Audio
- Dual Level Audio

Page 5 (P-5)

- Datalogging Mode
- Month
- Day
- Year
- Hour

- Minute
- Second

Note: Parameters on Page 5 (P-5) are not user editable using the instrument Setup Mode. To edit these parameters, please use the Ludlum LUMIC Calibration or datalogging Software.

Setup Parameter Adjustment: Tap (pressing the button and releasing under a half second) the MODE button to adjust the value for the selected item. When the appropriate value is selected for that item, tap the ON/ACK button to move to the next item. When the desired value is displayed, tap the UNITS button to advance to the next parameter. The graph will display the total number of parameters available on the current Setup Page, and the current parameter's position will be blinking. *When the Model 3001 is in PROTECT mode (dipswitch setting), the Setup parameters will cycle through to display the set values, but changes are not possible.*

The order of Setup parameters for the Model 3001 is as follows:

Setup Page 1

Calibration Constant Mantissa (Default 200) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Value is used with Calibration Constant Exponent to express counts per selected exposure unit. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (1-9)

Calibration Constant Exponent (Default 6) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Value is used with Calibration Constant Mantissa to express counts per selected exposure unit. Available values are:

Ones Place (0-9)

Tens Place (0-1)

Dead Time Correction (Default Varies) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Value is in microseconds. User should be able to zero when checking instrument with pulser or function generator. Otherwise, higher count rates will result in non-linearities. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

See the following table for examples with specific detectors.

Detector	Dead Time	Calibration Constant
Model 44-2	10 μ sec \pm 5 μ sec	100e+08
Model 44-6	100 μ sec \pm 35 μ sec	709e+05
Model 44-7	262 μ sec \pm 35 μ sec	139e+06
Model 44-9	85 μ sec \pm 35 μ sec	200e+06
Model 44-38	100 μ sec \pm 35 μ sec	709e+05

Dead Time Correction 2 Mantissa (Default 0) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (1-9)

Dead Time Correction 2 Exponent (Default 0) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Value between (-6 to -12)

Efficiency (Default 15.0%) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0 or 1)

Normally the efficiency is used on a per detector basis or 4pi basis, i.e. the efficiency of the detector is calculated by dividing the count rate received from a source by the total disintegration rate of the source. When either dpm or Bq units are chosen, the use of the 4pi efficiency allows the display of the source size or activity on the Model 3001 display.

But if the user desires to have the Model 3001 show results in terms of dpm/100 cm², the user could manipulate the efficiency to produce this result by multiplying the efficiency times the ratio of the detector area to 100 cm². For example, using a detector with an area of 15 cm², if we start with 10% efficiency to measure in dpm, then the parameter could be changed to 1.5% to measure in dpm/100 cm².

Or likewise for Bq/cm², efficiency could be calculated as: efficiency = count rate/disintegration rate*detector area (in cm²). For example, with the same detector as above with an area of 15 cm², and starting with an efficiency value of 15%, then the parameter could be changed to 225% to measure in Bq/cm².

High Voltage (Default 900 Volts) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Value is in Volts. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

k Multiplier (on/off)

Note: k multiplier also activates left-most decimal point. If k multiplier is used, Hundreds Place value is limited to 0 and 1.

Pulser Threshold (Default 35 millivolts) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Note: Audio is enabled for this parameter only for east of setting.

Detector Current Overload Threshold (Default 100) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Setting the Detector Current Overload Threshold to 0 disables the Current Overload Alarm. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

A long press (at least a half second) of the DETECTOR button will cause the display to show the actual live current reading. Repeating this process will return you back to the overload threshold setpoint.

Loss of Count Alarm Time (Default 60 seconds) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Setting the Loss of Count Alarm Time to 0 disables the alarm. Value is in seconds. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Setup Page 2

Primary RATE/MAX Units and Minimum Display (Default 0 μ R/hr) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Number of Decimal Places (0, 1, or 2)

Minimum Units (See list below.)

cpm	kcpm	cps	kcps	dpm	kdpm	Bq
kBq	μ R/h	mR/h	R/h	μ Sv/h	mSv/h	Sv/h

Primary Units RATE/MAX Maximum Display (Default 999 mR/hr) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Units will be the same as selected earlier with Primary Units. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (k on or off – cpm, cps, dpm, Bq; μ , m or none for R/h and Sv/h)

Primary Units RATE/MAX Mode Alarm Point 1 and 2 (Default 000) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Units will be the same as selected earlier with Primary Units. The ALARM LCD Segment will flash to indicate Alarm 1 and will be on solid to indicate Alarm 2. Set this Alarm to 000 to disable. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (k on or off – cps, cpm, Bq; μ , m or none for R/h and Sv/h)

Note: If the Primary Units has changed to a value other than that used to previously set this Alarm Point, the Alarm Point will be reset to 000.

Primary Count Units and Minimum Display (Default 0 μ R) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Number of Decimal Places (0, 1, or 2)

Minimum Units depend on the selected Primary Units (See list below.)

Primary Units	Primary Count Units Available	Primary Units	Primary Count Units Available
cps kcps	cps c	μ R/h mR/h R/h	μ R/h μ R mR/h mR R/h R
cpm kcpm	cpm c		
Bq kBq	Bq d	μ Sv/h mSv/h Sv/h	μ Sv/h μ Sv mSv/h mSv Sv/h Sv
dpm kdpm	dpm d		

Primary Count Alarm Point 1 and 2 (Default 000) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Units will be the same as selected earlier with Primary Count Units. Primary Count Units of c or d will not be displayed, but mR or μ Sv will. The ALARM LCD Segment will flash to indicate Alarm 1 and will be on solid to indicate Alarm 2. Set this Alarm Point to 000 to disable. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (k on or off – cpm, cps, dpm, Bq; μ , m or none for R/h and Sv/h)

Note: If the Primary Count Units has changed to a value other than that used to previously set this Alarm Point, the Alarm Point will be reset to 000.

Setup Page 3

Secondary RATE/MAX Units and Minimum Display (Default 0 cpm) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. Available values are:

Number of Decimal Places (0, 1, or 2)

Minimum Units (See list below.)

cpm	kcpm	cps	kcps	dpm	kdpm	Bq
kBq	μ R/h	mR/h	R/h	μ Sv/h	mSv/h	Sv/h

Secondary Units RATE/MAX Maximum Display (Default 999 kcpm) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. If the Secondary Units is off, this parameter will be skipped. Units will be the same as selected earlier with Secondary Units. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (k on or off – cpm, cps, dpm, Bq; μ , m or none for R/h and Sv/h)

Secondary Units RATE/MAX Mode Alarm Point 1 and 2 (Default 000) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. If the Secondary Units is off, this parameter will be skipped. Otherwise, units will be the same as selected earlier with Secondary Units. The ALARM LCD Segment will flash to indicate Alarm 1 and will be on solid to indicate Alarm 2. Set this Alarm Point to 000 to disable. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (k on or off – cpm, cps, dpm, Bq; μ , m or none for R/h and Sv/h)

Note: If the Secondary Units has changed to a value other than that used to previously set this Alarm Point, the Alarm Point will be reset to 000.

Secondary Count Units and Minimum Display (Default c) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. If the Secondary Units is off, this parameter will be skipped. Available values are:

Number of Decimal Places (0, 1, or 2)

Minimum Units depend on the selected Primary Units. (See list below.)

Secondary Units	Secondary Count Units Available	Secondary Units	Secondary Count Units Available
cps kcps	cps c off	μ R/h mR/h R/h	μ R/h μ R mR/h mR R/h R off
cpm kcpm	cpm c off		
Bq kBq	Bq d off	μ Sv/h mSv/h Sv/h	μ Sv/h μ Sv mSv/h mSv Sv/h Sv off
dpm kdpm	dpm d off		

Secondary Count Alarm Point 1 and 2 (Default 000) - Tap ON/ACK to select the value to adjust, and MODE to adjust the value. If the Secondary Units is off, this parameter will be skipped. Otherwise, units will be the same as selected earlier with Secondary Count Units. The ALARM LCD Segment will flash to indicate Alarm 1 and will be on solid to indicate Alarm 2. Set this Alarm to 000 to disable. Available values are:

Ones Place (0-9)

Tens Place (0-9)

Hundreds Place (0-9)

Number of Decimal Places (0, 1, or 2)

Range (**k** on or off – cpm, cps, dpm, Bq; μ , m or none for R/h and Sv/h)

Note: If the Secondary Count Units has changed to a value other than that used to previously set this Alarm Point, the Alarm Point will be reset to 000.

Setup Page 4

Response Time (Default 0 - auto) – Tap ON/ACK to select the value to be adjusted and MODE to adjust the value. Setting the Response Time to a fixed value is useful primarily when performing surveys to a fixed MDA (Minimum Detectable Activity) level. Setting the Response Time to 0 will enable the Auto-Response mode for the Model 3001 (see the next parameter). Available values for the Response Time (in seconds) are:

Ones Place (0-9)

Tens Place (0-6, 6 forces max response time of 60)

Auto-Response Rate (Default S) - Tap MODE to select Fast (**F**) or Slow (**S**). When operating in Auto-Response mode, the Model 3001 will vary the Response Time based on the Auto-Response Rate selected (**Fast** or **Slow**) and the current Count Rate. The following table shows the response time for different count rates when these auto response modes are chosen:

Count Rate	Auto Response Time – Fast (Seconds)	Auto Response Time – Slow (Seconds)
Less than 3 kcpm (50 cps)	10.5	21
Between 3 kcpm and 4 kcpm (67 cps)	8.4	16.8
Between 4 kcpm and 6 kcpm (100 cps)	6.3	12.6
Between 6 kcpm and 12 kcpm (200 cps)	4.2	8.4
More than 12 kcpm	2.1	4.2

The Model 3001 also utilizes a Step function in Auto Response mode, which enables faster response to a significant increase or decrease in Count Rate. When the instrument detects a sudden change in count rate from the detector, the response time is reduced to 1 second to quickly show the new value.

Operational Modes (Default 0 - All modes available) - Tap MODE to adjust the value. Available values are:

- 0 – RATE, MAX, and COUNT Modes
- 1 – RATE and MAX Modes only
- 2 – RATE and COUNT Modes only
- 3 – RATE Mode only

Count Time (Default 1 minute) - Tap ON/ACK to select the value to adjust and MODE to adjust the value. Setting Count Time to 0 enables continuous count until reset. If 9 minutes are selected, then the maximum seconds value is 60; otherwise, the maximum seconds value is 59. Available values are:

Ones Place (0-9)

Tens Place (0-6, 6 only available if minutes value is 19)

Hundreds Place (0-9)

Note:

The UNITS button can be used to advance to the next parameter. To send Setup mode and save the current setting, press and hold the UNITS button for approximately 5 seconds, or press the UNITS button for at least a half a second two times.

Auto Shutdown Time (Default 0 - off) - Tap MODE to adjust the value. Setting Auto Shutdown Time to 0 disables Auto Shutdown. Available values are:

Ones Place (0-9) hours

Backlight Threshold (Default LO –most sensitive) - Tap MODE to adjust the value. Available values are:

LO – Backlight comes on at a higher ambient light level

There are two different options for **Backlight Threshold LOW**, which can be set through software. Locate the Parameters in the Others tab. Enter the desired number into the “Value” field.

- **Low Light Turn On:** The value will set the threshold that the backlight will turn on at. (This should be lower than the Turn Off value.)
- **Low Light Turn Off:** The value will set the threshold that the backlight will turn off at. (This should be higher than the Turn On value.)

HI – Backlight comes on at a lower ambient light level

There are two different options for **Backlight Threshold HIGH**, which can be set through software. Locate the Parameters in the Others tab. Enter the desired number into the “Value” field.

- **High Light Turn On:** The value will set the threshold that the backlight will turn on at. (This should be lower than the Turn Off value.)
- **High Light Turn Off:** The value will set the threshold that the backlight will turn off at. (This should be higher than the Turn On value.)

OFF – Disables backlight

Sigma Audio (Default OFF) - Tap MODE to adjust the value. Available values are:

ON – Sigma Audio Enabled

OFF – Sigma Audio Disabled (Normal 'Click' Audio active)

Note

The Sigma Mode may not work properly until the 8-second background is completed. If the 8-second background count is too low or too high, the device will continually beep until an acceptable background count rate (500 cpm - 200 kcpm) is obtained. Sigma Audio is disabled in Count Mode. If Count audio and Sigma Audio are enabled, Count audio will be "click" audio.

Dual Level Audio (Default OFF) - Use MODE to adjust the value. Available values are:

ON – Dual Level Audio Enabled (Audio settings: Hi, LO, and Mute)

OFF – Dual Level Audio Disabled (Audio settings: Hi and Mute)

Setup Page 5

Note:

Parameters on Page 5 (P-5) are not user editable using the instrument Setup Mode. To edit these parameters, please use the LMI Lumatic Datalog Software.

Datalogging Mode (Default 1) – Displays the current Datalogging Mode selection for the device. Valid values are: 1-3

Month (Default 1) – Displays the month of the Real Time Clock. Available values are: 1-12

Day (Default 7) – Displays the numerical day of the Real Time Clock. Available values are: 1-31

Year (Default 14) – Displays the decade and year of the Real Time Clock. Available values are: 00-99

Note: Valid year range is 2000-2099.

Hour (Default 16) – Displays the hour of the Real Time Clock in 24-hour format. Available values are: 0-23

Minute (Default 30) – Displays the minute of the Real Time Clock. Available values are: 00-59

Second (Default 00) – Displays the second of the Real Time Clock. Available values are: 00-59

Section**5**

Datalogging

The datalogging feature of the Model 3001 allows the user to log radiation readings with the use of a handle-mounted LOG button. Data can be logged in any of the Model 3001's operational modes (RATE, MAX, and COUNT). Up to 1000 data points can be taken and stored internally. Depending on the chosen Datalogging Mode, the user can quickly save logged data using a single Location ID, or select from up to 250 stored Location IDs. A Real Time Clock (RTC) is utilized to time and date stamp each datalog entry.

The saved log data and stored location IDs can be viewed, downloaded, and erased using the LMI Lumic software kit, which includes a USB cable. Editing of datalogging parameters and RTC via setup mode on the instrument is disabled.

Setup of the datalogging parameters is also done through the LMI Lumic software. The advanced user or administrator can set:

- Datalogging Mode
- Real Time Clock (RTC) Values
- Location IDs

Datalogging Operation – Mode 1

Datalogging Mode 1 will store the logged data using only the first Location ID in the Location ID table.

- When the LOG button is pressed, the current radiation reading and other log data is saved along with the first Location ID in the format specified in the Specifications at the end of this section.
- In COUNT Mode:
 - With a set count time, the LOG button is disabled until completion of a Scaler count.
 - For a continuous Scaler count (Scaler time is 0), the LOG button is enabled at all times.
- For approximately 2 seconds, the LCD will display the Datalog Table index for the newly saved log data.



Figure 12: Display showing a Datalog Table index of 12.

- After displaying the Datalog Table index, the instrument will return to the previous mode of operation.
- The Datalog Table is linear; once the table is full, no further writes will be allowed until the Datalog Table is erased.
- If an attempt is made to write to a full Datalog Table, a beep will sound and the maximum Datalog Table index (1.00 k) will be displayed for about 2 seconds. The instrument will then return to the previous mode of operation.

Datalogging Operation – Mode 2

Datalogging Mode 2 will allow the user to choose the Location ID (by Location ID Table index) to store with the logged data.

- When the LOG button is pressed, the current radiation readings and other log data are temporarily stored.
- In COUNT Mode: (Datalogging Mode 1 is always utilized in COUNT Mode.)
 - With a set count time, the LOG button is disabled until completion of a Scaler count.
 - For a continuous Scaler count (Scaler time is 0), the LOG button is enabled at all times.
- The LCD display will show a possible Location ID Table index for the user. The index will be auto-incremented from the previously used index.



Figure 13: Display showing a Location ID Table index of 36.

- The ones digit of the index will be blinking, indicating a changeable value. The user may then enter the preferred Location ID Table index by tapping the ON/OFF (to select the digit) and MODE (to change the value) buttons as in Setup mode.
- Once the user has the preferred Location ID Table index entered, tap the LOG button to save the log data.
- For approximately 2 seconds, the LCD will display the Datalog Table index for the newly saved log data. (See Figure 12.)
- After displaying the Datalog Table index, the instrument will return to the previous mode of operation.
- The Datalog Table is linear; once the table is full, no further writes will be allowed until the Datalog Table is erased.

Datalogging Operation – Mode 3

Data Logging Mode 3 will automatically record data log records using the current location with the user settable auto log interval (settable to record every 1 to 1800 seconds).

- When the Device is turned on and in RATE or MAX mode, the current radiation reading and other log data is recorded at the set interval.
- In SCALER Mode:
 - Auto log will be paused and will only log at the completion of each count.
 - For a continuous Scaler count (Scaler time is 0), logs will be recorded at auto log interval.
- The Data Log Table is linear; once the table is full, no further writes will be allowed until the Data Log Table is erased.

Note:

Tapping the LOG button will start a count and automatically log the result. Long-pressing the LOG button during the count will toggle the display between the timer and reading.

Specifications

Handle mounted LOG button

Real Time Clock (RTC)

Internal Storage of up to 250 32-byte Location IDs

Internal Storage of up to 1000 64-byte Datalog records

Datalog format:

- Format Version (1 Byte)
- Month (1 Byte)
- Year (2 Bytes)
- Day (1 Byte)
- Hour (1 Byte)
- Minutes (1 Byte)
- Seconds (1 Byte)
- Logged Value (4 Bytes)
- Range (1 Byte)
- Units (1 Byte)
- Mode (1 Byte)

- Detector Number (1 Byte)
- Status (1 Byte)
- Reserved (2 Bytes)
- Elapsed Count Time in seconds (4 Bytes)
- Scaler Count Time in seconds (4 Bytes)
- Location ID (32 Bytes)

Section
6**Software**

Connecting to Software

The Model 3001 software are sent with a standard 2-meter cable (A 5-meter cable can be provided if requested. However, any cable longer than 2 meters may have issues with some USB-hubs/computers.)

To connect an instrument to the computer, please connect one end of the USB cable to the instrument first, then the other end to the computer. Do not connect both ends to the computer.

Please allow Windows® a moment to install the proper HID drivers for the instrument before trying to use any software.

Note:

We recommend that you plug the USB cable into the back of your PC that connects to your motherboard instead of a USB hub.

Note:

Some parameters may only be edited in software, such as the backlight thresholds, COUNT Display Mode, COUNT Audio Mode, and Setup prtct.

Extended Features

This section describes features that have been added to the firmware since the last software release. If the version of your Lumic Calibration Software is higher than 1.2.4.4, these parameters may have been moved to the Device tab.

Device Default Detector

- -1 – Last used detector
- 0 – Boot up using Detector 1
- 1 – Boot up using Detector 2
- 2 – Boot up using Detector 3
- 3 – Boot up using Detector 4

Detector 1 Channel 1 CPSOffset


- Counts per second to be subtracted from base counts.
- Applied before any deadtime correction or unit conversion.

Device Button Handle (operating mode) (button press type)

- The handle button can functionality can be mapped to do different things. There are 6 variables to configure this, 3 types of button presses for each of the 2 types of operating modes (listed below).
 - RateMap – Applies when operating in Rate or Max Rate Modes.
 - CntMap – Applies when operating in Count Mode.
- 1 – 3 defines different button press types.
 - 1 – Press and release in less than 0.5 seconds (tap)
 - 2 – Press and release 0.5 – 2 seconds (long press)
 - 3 – Press and hold for > 2 seconds (very long press)
- The value of the parameter defines the action taken. The available actions are as follows:
 - 0 – Do nothing
 - 1 – Take Log
 - 2 – Change units (like short press of Units button)
 - 3 – Toggle click audio
 - 4 – Change Modes (like short tap of Mode button)
 - 5 – Unused

- 6 – Toggle display of Count Time and Value while counting in count mode.
- 7 – Reset rate while in rate mode
- 8 – Advance count mode state. Starts count when ready, resets counter while running or when finished.
- 9 – Cycles detectors (like short tap of detector button. 3001 only)
- Any other number will be the same as 0
- As an example, the settings in the picture below will produce the following actions if the Log (handle) button is held for the listed amount of time.
 - If in rate mode:
 - < 0.5 sec – Change units (note no log will be taken)
 - 0.5 – 2 sec – Reset Rate to 0
 - 2 sec – Take a log
 - If in count mode:
 - < 0.5 sec – Change units (note no log will be taken)
 - 0.5 – 2 sec – Toggle the display of count value and count time while a count is running.
 - 2 sec – Count mode will start counting if currently reset, or will reset if currently counting or counting is finished.

Device Button Handle RateMap 1	2	<input type="checkbox"/>
Device Button Handle RateMap 2	7	<input type="checkbox"/>
Device Button Handle RateMap 3	1	<input type="checkbox"/>
Device Button Handle CntMap 1	2	<input type="checkbox"/>
Device Button Handle CntMap 2	6	<input type="checkbox"/>
Device Button Handle CntMap 3	8	<input type="checkbox"/>

A dark gray square containing the word "Section" in a bold, sans-serif font at the top, and a large, white, bold number "7" in the center.

Advanced Features

Dead Time Correction

All pulse counting detectors have a “dead time” in which the detector is unable to register another event. In relatively low fields this is not an issue. However, as the field strength approaches the high end of the detector’s range, dead time causes the pulse rate to become non-linear with respect to the real radiation field. Dead time correction is used to linearize the measurements, allowing a wider linear response range for a given detector.

This instrument uses a second order dead time correction using the equation

$$\text{rate}_c = \frac{\text{rate}_m}{[1 - (DTC_1 * \text{rate}_m) + (DTC_2 * \text{rate}_m^2)]}$$

where DTC_1 and DTC_2 are the dead time correction coefficients, rate_m is the measured count rate, and rate_c is the corrected rate.

Dead time correction coefficients are configured through the device setup menu or through *Lumic Calibration* software. *Lumic Calibration* software also includes a wizard that will automate finding and setting the correct coefficients.

DTC_1 is represented in micro seconds (μs) on the device setup menu and in the software. The equation must be calculated in terms of seconds (s) and must be scaled appropriately. DTC_2 has more complex units and are not shown. However, the representation of the coefficient on the setup menu and in software do not require scaling.

Note:

Setting the DTC_2 to zero will disable the second order correction, while setting DTC_1 to zero will disable both, regardless of the value of DTC_2 .

Units

Depending on the chosen display units, different features will affect value of the reading. The following table lists the features that apply to each of the display units.

Unit	Feature
cpm cps counts	Dead time correct
Bq dpm disintegrations	Dead time correction Efficiency
R/h R	Dead time correction Calibration constant
Sv/h Sv	Dead time correction Calibration constant R to Sv conversion

These apply to all modes.

R to Sv Conversion

The **R to Sv conversion** is a setting available in *Lumic Calibration* software. It defines the conversion factor between R and Sv. Since the **calibration constant** is directly tied to R (counts per R), calibrating Sv requires a correct setting of both **calibration constant** and **R to Sv conversion**.

Example: A **R to Sv conversion** factor of 0.0106 will cause a reading of 10.6 mSv/h while in a 1 R/h field.

Software Calibration Tools

Lumic Calibration software includes wizards that will assist in calibrating and plateauing detectors. After configuring the wizard for a specific detector, the wizard will automate much of the data collection and calculation required for calibration.

Instrument Calibration

This instrument provides the ability to digitally enter a voltage for both high voltage and threshold. As such, the instrument itself needs to be calibrated. This happens at the factory when the instrument is built and should not require re-calibration. However, if calibration is required (due to board rework, etc.) the calibration settings are available in *Lumic Calibration* software.

Other Device Data

The following parameters on the instrument allow recording import device information within the device:

Firmware Version: This is a read-only presentation of the firmware version. With a version of 4983xnny.zzzz, the xyy will show up on the device screen during the power-on sequence and signifies the released version.

Device – Model Name: This should match the model name on the front face of the instrument.

Device – Serial Number: This should match the serial number of the instrument.

Detector – Model: This can store the model of the detector the instrument was calibrated for.


Detector – Serial Number: This can store the serial number of the detector the instrument was calibrated for.

Real-time Streaming

Lumic Calibration software and this instrument have the ability to stream data from the instrument to a computer. The data can be viewed live inside software or can be recorded on file.

Multiple user-selected parameters can be streamed simultaneously including:

- Remote display of the screen
- All unit values
- Device status
- Live HV current measurement

A dark gray square with the word "Section" in white at the top and a large white number "8" in the center.

Safety Considerations

Environmental Conditions for Normal Use

Indoor or outdoor use (While rain resistant, user is cautioned to avoid getting water through detector opening.)

No maximum altitude

Temperature range of -20 to 50 °C (-5 to 122 °F) , may be certified for operation from -40 to 65 °C (-40 to 150 °F)

Maximum relative humidity of less than 95% (non-condensing)

Pollution Degree 3 (as defined by IEC 664): (Occurs when conductive pollution or dry nonconductive pollution becomes conductive due to condensation. This is typical of industrial or construction sites.)

Not certified for use in an explosive atmosphere

Warning Markings and Symbols

The Model **3001** Survey meter is marked with the following symbols:



CAUTION (per ISO 3864, No. B.3.1): designates hazardous live voltage and risk of electric shock. During normal use, internal components are hazardous live. This instrument must be isolated or disconnected from the hazardous live voltage before accessing the internal components. This symbol appears on the side panel. Be sure to take the precautions noted in the next section whenever necessary.



The “**crossed-out wheeie bin**” symbol notifies the consumer that the product is not to be mixed with unsorted municipal waste when discarding. Each material must be separated. The symbol is placed on the label located on the side panel. See section 7, “Recycling,” for further information.



The “CE” mark is used to identify this instrument as being acceptable for use within the European Union.


Caution!

The operator or responsible body is cautioned that the protection provided by the equipment may be impaired if the equipment is used in a manner not specified by Ludlum Measurements, Inc.

Cleaning and Maintenance Precautions

The Model 3001 may be cleaned externally with a damp cloth, using only water as the wetting agent. Observe the following precautions when cleaning or performing maintenance on the instrument:

1. Turn the instrument OFF and remove the batteries.
2. Allow the instrument to sit for one minute before cleaning the exterior or accessing any internal components for maintenance.

A dark gray square with the word "Section" in white at the top and a large white number "9" in the center.

Revision History

NOTE: This section of the manual will be updated with each revision of the Model 3001 in order to document changes over time. Ludlum Measurements' policy is to provide free software upgrades to instruments for the life of the instrument.

July 2015: New manual.

November 2016: Added description of headphone option. Corrected the location of the low-battery indicator. Updated the table on page 4-13. Removed "Units selected as Primary will not be available as Secondary" on page 4-11. Updated Sigma Audio Note on page 4-14. Added information about rate reset on page 2-6. Added different options for protect mode to firmware on page 4-3. Took out references to Triad software and replaced with Lumic. Added information about the new parameters Backlight Threshold, Device COUNT DisplayMode, Device COUNT AudioMode, Dead Time Correction 2. Added Section 5 Connection to Software, to most PCs.

April 2017: Corrected references of Audio button to Detector button. Firmware was updated, which affects Section 2 with more screens at turn-on and changes in what buttons do and how to press/tap them. Updated Audio in Section 3 Specifications. In Section 4 updated Default Values including adding P5 values, added a note on page 4-3 before Setup Protect explanation, added sentence about backlight configuration in note at end of Entering Setup Mode section on page 4-4, detailed types of presses/taps in setting up pages, changed values of Dead Time Correction 2 Exponent to -6 to -12, added Dual Level Audio default at end of Setup page 4, updated Response times in table on page 4-16. Added Section 5 Datalogging. Section 6 updated to Software with Connecting to Software and Extended Features as sub-sections. Added Section 10 Options. Added Appendix A Typical Values. Added front-panel drawing 498 x 409M in front of manual.

November 2018: Datalogging is no longer an "option" but a standard features with the logging button as standard on the handle, updated the photo on inside cover, added a paragraph at the end of Section 1, and deleted 3000 Series Datalogging Option in Section 11. Updated IP Rating to 33 from 65, based on further testing, Updated Section 4, table now Setup Parameters and

explanations of pages updated, added P5, updated alarm option to two alarms. Made corrections to first paragraph of Section 6, under Connecting to Software. Added Advanced Features as Section 7. Appendix A now Auxiliary Communications. Updated Drawing 498 x 409M. Added serial # and firmware # to cover.

Section 10 Recycling

Ludlum Measurements, Inc. supports the recycling of the electronics products it produces for the purpose of protecting the environment and to comply with all regional, national, and international agencies that promote economically and environmentally sustainable recycling systems. To this end, Ludlum Measurements, Inc. strives to supply the consumer of its goods with information regarding reuse and recycling of the many different types of materials used in its products. With many different agencies – public and private – involved in this pursuit, it becomes evident that a myriad of methods can be used in the process of recycling. Therefore, Ludlum Measurements, Inc. does not suggest one particular method over another, but simply desires to inform its consumers of the range of recyclable materials present in its products, so that the user will have flexibility in following all local and federal laws.

The following types of recyclable materials are present in Ludlum Measurements, Inc. electronics products, and should be recycled separately. The list is not all-inclusive, nor does it suggest that all materials are present in each piece of equipment:

Batteries	Glass	Aluminum and Stainless Steel
Circuit Boards	Plastics	Liquid Crystal Display (LCD)

Ludlum Measurements, Inc. products that have been placed on the market after August 13, 2005, have been labeled with a symbol recognized internationally as the “crossed-out wheelie bin,” which notifies the consumer that the product is not to be mixed with unsorted municipal waste when discarding. Each material must be separated. On the Model 3001, the symbol will be placed on the serial number label located on the side of the instrument.

The symbol appears as such:



Section
11**Options**

Lumic Calibration Kit (part # 4498-1018): The kit includes calibration software plus the cable required for calibration. The software allows users to collect data and read, print, and save device parameters. It allows administrators to adjust device parameters from one device to another.

Lumic Datalogging Kit (part # 4498-1019): The kit includes datalogging software plus the required cable. The software allows users to collect data and read, print, and save device parameters. It allows administrators to adjust device parameters from one device to another.

Headphone Option (part # 4498-555): This provides the Model 3000 series of instruments with a jack and circuitry required for a standard headphone plug. Ludlum Measurements also offers mono/stereo headphones with volume control.

Shoulder Strap (part # 4498-868): This adjustable, padded strap comes with a kit to attach it to a Model 3000 series instrument.

Protective Storage/Transport Case (part # 2312958): This is a medium-sized, foam-padded, rugged case that can be secured with a padlock. It is fitted with a manual pressure relieve valve for air transport, providing water and dust-proof protection for sensitive instruments.

Bluetooth 4.0 LE® (*Bluetooth Low Energy, sometimes referred to as Bluetooth Smart*) (part # 4498-1024): This feature permits wireless transmission of readings from a connected instrument, allowing operators to remotely monitor the live data on the screen of their mobile device. When paired with Ludlum's Lumic Linker App, the operator can also seamlessly send data to the *RadResponder Network, which provides a central location for up-to-date information from operators in the field. Reported data includes user, radiometric survey, survey notes, and GPS location, as well as details about the instrument and detector being used. (See the Model 3000 Series Bluetooth LE® Addendum for more details on how to use this option.)

Section
12**Standard Parts List**

Part Description**Part #**

Model 3001 Digital Survey Meter
Model 3001 Main Board
LCD 82 mm x 61.64 mm
Unimorph with Harness
Model 3001 Bezel Assembly
Main Keypad Membrane Switch
Battery Holder Modified

48-4036
5498-440
7498-417
8498-542
4498-476
7498-425
7498-458

Optional

Model 3001 Headphone Option
Lumic Calibration Software
Lumic Datalog Software

4498-555
4498-578
4498-775



Auxiliary Communications

AuxComm Overview

AuxComm, short for Auxiliary Communications, is a feature included on certain Ludlum instruments. An AuxComm port allows the instrument to expand its capabilities with a variety of external devices through a standard serial interface.

Requirements

In order to take advantage of the AuxComm functionality, you will need the following:

- An instrument with Printed Circuit Board Number (PCB#)
PCB# 5519-042 (Model 3000, 3001)
PCB# 5519-063 (Model 3002)
PCB# 5498-770 (Model 3276 **375 Modes only**)
- An AuxComm enabled firmware
- A wiring harness to connect the target device
- The correct mode set for the AuxComm port

Note:

Checking PCB: If there is a need to verify the correct PCB, open the instrument to look at the number printed on the PCB. (See Figure 14.)




Figure 14: Example – Printed Circuit Board (PCB) number location.

Usage

Usage of the AuxComm port is generally very simple from the user's perspective. Most modes only require everything to be plugged in and then should “just work.” However, there are a couple of caveats.

If using the SLURM protocol, the user will need to pair the instrument with a mobile device using the hot-key listed below.

AuxComm Usage - Hot-key

The hot-key is triggered by simultaneously pressing the UNITS and MODE buttons as seen in Figure 15.

There are two types of hot-key presses:

- 1 second:** If the hot-key combination is held for a little over 1 second (but less than 2), the AuxComm port will toggle on and off each time this is done.

>2 second: If the hot-key combination is held for longer than 2 seconds, the instrument will present the user with a PIN needed for pairing. See Figure 16.



Figure 15: Example – AuxComm Hot-key.

Caution!

The pairing PIN is only used for SLURM.

AuxComm Usage - Pairing

“Pairing” is the process to connect an instrument with a mobile device. It helps verify you are connecting to the correct devices, as well as configures the encryption. The pairing process should only be done when both the instrument and mobile device are ready. Ensure the App or software is installed, and then follow these steps:

1. Open connection window in App/Software.
2. Use 2-second PIN hot-key to get PIN. This also ensures the wireless module is turned on.
3. Select “Scan” in App/Software.
4. When instrument is found by mobile device, select device in list and then “Pair.”
5. Enter PIN in App/Software from step 2.




Figure 16: Example – Device Selection In Lumic Linker App (iOS)




Figure 17: Example – Model 3001 Pairing




Figure 18: Example – Model 3001 Found in Lumic Linker App (iOS)




Figure 19: Example – Entering PIN in Lumic Linker App (iOS)

Settings

The settings for AuxComm are only configurable through software. See Figure 20.




Figure 20: Example – Lumic Calibration Software AuxComm Menu

AuxComm Setting - Enabled

This setting enables or disables the AuxComm port. If this is set to “Off,” the port is completely disabled and will stay that way until set to “On.”

Note:

When AuxComm is ”Off” the hot-key to turn AuxComm on and off on the fly will not function.

Table A-1: AuxComm Enabled

Setting	State
0	Off
1	On

AuxComm Setting - Mode

The mode setting for AuxComm configures how the port will function. Table A-2 lists the modes and their baud rates.

Table A-2: AuxComm Mode

Setting	Mode	Baud Rate
0	SLURM	115200
1	LMI Direct	115200
2	375 Standard	2400
3	375 Ethernet	2400
4	375 for 272D	2400
5	Simple Counts	115200

SLURM: Used to encapsulate “LMI Direct” messages with encryption and checksums. This mode is primarily used for Bluetooth®.

LMI Direct: Full exposure to the “LMI Comm” Communication protocol. This allows reading and writing all parameters, data streaming, etc.

375 Standard: Legacy support for Model 375 Standard output. This is compatible with the Model 375 Logger software and some versions of the Model 3276 Ethernet Kit (LMI# 4498-958).

375 Ethernet: Extended version of “375 Standard” that is intended for use with Model 3276 Ethernet Kit (LMI# 4498-958 **FW# N09 Only**) and the Ludlum Universal Software.

375 for 272D: Special version of “375 Standard” that makes the 272D display follow the instrument display more closely with regards to leading zeros and decimal places.

Simple Counts: Mode that automatically transmits raw, uncorrected counts. The count interval is based on the AuxComm Auto Mode Interval variable. The output value is the cumulative counts that occurred in the previous interval.

Note:

The Model 3276 is the only unit that is wired to properly interface with the 375 Ethernet module.

AuxComm Setting - Auto Off Time

This setting defines the number of minutes the AuxComm port will sit idle before turning off. The primary use for this parameter is to save power while using Bluetooth® by turning off the Bluetooth® module if there is no connection. Table A-3 shows valid settings.

Note:

The Model 3276 is the only unit that is wired to properly interface with the 375 Ethernet module.

Table A-3: AuxComm Auto Off Time

Setting	Description
0	Disabled
1-9	Inactivity minutes before auto off

AuxComm Setting - Power on at Boot-up

This setting determines whether or not the AuxComm module is turned on when the instrument is powered on. In addition to power savings, this setting is also intended for use when a user has a need for “radio silence” (no RF broadcasts). This setting does not override the hot-key (Figure 15) function. Again, the primary use for this parameter is while using Bluetooth®. Table A-4 shows valid settings.

Note:

This setting should be 0 when using any of the 375 modes.

Table A-4: AuxComm Power on at Boot-up

Setting	Description
0	AuxComm On at boot-up
1	AuxComm Off at boot-up

AuxComm Setting - Write Protect

This setting is an additional security measure that prevents commands that would modify parameters from being executed over the AuxComm port. While using the instrument equipped with a wireless radio, such as the Bluetooth® module, there is the potential for the instrument becoming mis-configured, inadvertently or intentionally, from a remote location. This setting is meant to ensure that cannot happen. Table A-5 shows valid settings.

Note:

This has no effect in 375 modes.

Table A-5: AuxComm Write Protect

Setting	Description
0	AuxComm Write Protect Disabled
1	AuxComm Write Protect Enabled

Caution!

If write protect is enabled, it can only be disabled by physically connecting it to the instrument through the USB port.

AuxComm Setting - Encryption

This setting controls whether or not encryption is enabled on the AuxComm port. While using the instrument equipped with a wireless radio, such as the Bluetooth® module, there is the potential for external eavesdropping or attack from a third party. By utilizing an industry standard encryption, this becomes much more difficult. Table A-6 shows valid settings.

Note:

This setting only applies to SLURM mode.

Table A-6: AuxComm Encryption

Setting	Description
0	AuxComm Encryption IS NOT used
1	AuxComm Encryption IS used

Note:

SLURM mode will still require the “pairing” sequence to connect to the instrument.

AuxComm Setting - Auto Mode Interval

This setting controls the number of seconds for each transmitted value in automatic stream modes. Table A-7 shows valid settings.

Note:

This setting only applies to Simple Counts mode as of this writing.

Table A-7: AuxComm Auto Mode Interval

Setting	Description
0	As fast as possible (1/4 second on Model 3000 family)
1-255	n-Seconds

Settings on Particular Models

The following table, Table A-8, shows the list of AuxComm settings available on various models.

Model	AuxComm					
	0	1	2	3	4	5
3000	x	x	x	x	x	x
3001	x	x	x	x	x	x
3002						
3005						
3019						
3276			x	x	x	