

Thermo Fisher Scientific
Personal DataRAM Communication Software version 2.10

OVERVIEW

The personal DataRAM communication software is for use with the pDR1000 and the pDR1200. Version 2.10 provides 64-bit operating system compatibility, but does not include graphing capability.

NOTE For older Operating Systems (pre 64-bit operating systems), you can use version v1.61 for its graphing capability. Please go to the pDR COM v1.61 directory for INSTALL.EXE file to load this software.

A. Install v2.10 pDRCOM Software

This installation procedure is described below:

1. Run pDRCOM0210setup.exe

RELEASE NOTES

Version 2.10 changes relative to version 1.61:

1. The graph data feature is no longer available.
2. .pdr data file format changed to human readable text. Data files extracted using the older version of pDR-Com will not open in this version of the program.
3. Improper spaces in data header were removed.